

Advanced, multi-user blood pressure monitor with enhanced features that provide trusted accuracy and reliable measurements every time.

Product Benefits

- Experience accurate readings with advanced cuff fit detection, avoiding errors from incorrect cuff placement
- Ensure precise measurements with integrated movement detection, minimizing the risk of inaccurate results due to movement
- Displays the average of all measurements providing a quick snapshot
- Easily track blood pressure over time with 60 measurement memory per user eliminating the need for manual recording

Comfort and Convenience

- Advanced inflation and deflation technology provides faster and more comfortable readings
- Up to 4 users provides convenience, cost savings, and enhanced data management capabilities
- Includes a carrying case and batteries for immediate use and portability capabilities

The A&D Promise

- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value



60 Measurement Memory

Irregular heartbeat symbol (IHB)/ AFiB

Cuff Fit Error

Movement Error

Pressure Rating Indicator

Battery Life Indicator



Clock Display

User



Hypertension CANADA

SET Button

Start Button



Unlock the Potential
Download the FREE
A&D Heart Track
Companion App
for manual entry



Technical Specifications

Measurement Method	Oscillometric
Memory	60 measurements for four users
Cuff Size	22-42 cm (8.6-16.5")
Measurement Range	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats / minute
Measurement Accuracy	Pressure: ± 3 mmHg, Pulse: $\pm 5\%$
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	14 x 10 x 6 cm (5.5 x 4.1 x 2.3")
Monitor Weight	255 g (9.0 oz) without batteries
Display Size	4.3 x 4.8 cm (1.7 x 1.9")
Clinical Test	ISO81060-2: 2013
Operating Conditions	+10°C to +40°C (50°F to 104°F) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-20°C to +60°C (-4°F to 140°F) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-767FAM Individual Unit (Box)	13 x 9.9 x 20 cm (5.1 x 3.9 x 7.7")	0.8 kg (1.7 lb)	0.0014 cm ³ (0.05 ft ³)	0-93764-60347-1
UA-767FAM Case (10 units)	51 x 28 x 23 cm (20 x 11 x 9")	8.0 kg (17.6 lb)	0.031 cm ³ (1.08 ft ³)	1-00-93764-60347-8

Optional Accessories

Name	Order Code
AC Power Adapter	TB-233

Compatible Replacement Cuffs- Packaged in individual box

Medium Cuff 23-37 cm (9.0-14.6")	UA-290A
Large Cuff 31-45 cm (12.2-17.7")	UA-291A
Wide Range Cuff 22-42 cm (8.6-16.5")	UA-420A
Cuffs are not made with natural rubber latex	

Wholesaler

	Codes
McKesson	082964
Kohl & Frisch	152971
LPG	35251
Imperial	201767

Included in the Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff 22-42 cm (8.6- 16.5")
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 1 AC adapter

Health Canada license - 94753



6695 Millcreek Drive, Unit 6
Mississauga, ON L5N 5R8, Canada
1-800-461-0991
andmedical.com