




Starting July 1st, Medi-Cal is adding a wider selection of A&D Medical blood pressure monitors and cuffs to meet the diverse needs of the Medi-Cal population.

With A&D Medical, your Medi-Cal beneficiaries can focus on their health journey with easy to use devices and the highly rated FREE A&D Heart Track app to log and trend their progress.



Blood Pressure Monitors Built to Last

- Simple to operate
- Extended 5 year warranty on the device
- Clinically validated with proven accuracy
- Broad selection to meet the needs of your diverse Medi-Cal beneficiaries
- Devices compatible with additional cuffs for a variety of sizes and styles

Best in Class Heart Track App




- Advanced analytics
- Integrated patient education
- Highly rated and easy to use
- Quickly log and share data
- Available for download by all users
- Manual entry available for non-wireless devices



FREE
Download



Let your Medi-Cal Beneficiaries select a customized blood pressure based on their needs. Monitors always include **5 year warranty**, **clinical validation**, and **one button operation**.

	 One Touch UA-660	 Multi User UA-770BLE	 Voice Guided UA-1040TBLE
Included Cuff Type (cm/in)	Wide Range 8.6-16.5" / 22-42 cm	Wide Range 8.6-16.5" / 22-42 cm	Wide Range 8.6-16.5" / 22-42 cm
Compatible with Different Cuff(s)	✓	✓	✓
Reading Memory	60	60/user	90
Averaging	✓	✓	Advanced
Irregular Heartbeat Detection	✓	✓	✓
Body Movement and Cuff Fit Error		✓	✓
Blood Pressure Rating Indicator	✓	✓	✓
Bluetooth® Connectivity		✓	✓
What is Included	4 AA Carrying Bag	4 AA Carrying Case	4 AA Carrying Case
Time & Date Recording	✓	✓	✓
Multiple Users		4 Users	
Verbal Assist			✓
MSRP	\$39.99	\$59.99	\$79.99
Recommended MAC	\$35.99	\$53.99	\$71.99
MediCal Pricing	\$20.80	\$34.94	\$45.00

**Listed blood pressure monitors without integrated cuffs can be operated using batteries (AA or AAA) or an AC Adapter (AC).

To place an order:

- E-mail custsat@andonline.com, fax 408-635-2338 or call 888-726-4772.
- Refer to account number 105026 when ordering.
- Orders require prepayment unless there is an approved credit application.

