

A premium Bluetooth® blood pressure monitor with advanced features and voice guidance, making it ideal for a wide range of users, including those who are visually impaired.

Product Benefits

- Audible instructions and results in English, Spanish and French for enhanced accessibility
- TriCheck™ minimizes the impact of individual measurement fluctuations by automating 3 sequential measurements and averaging in accordance with clinical guidelines
- Experience accurate readings with advanced cuff fit detection, avoiding errors from incorrect cuff placement
- Ensure precise measurements with integrated movement detection, minimizing the risk of inaccurate results due to movement
- Automatically records and tracks measurements with the highly rated Heart Track app to enable trending and sharing for better health management

Comfort and Convenience


- Advanced inflation and deflation technology provides faster and more comfortable readings
- Includes a carrying case and batteries for immediate use and portability

The A&D Promise

- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value



- Advanced Averaging
- 90 Measurement Memory
- Irregular Heartbeat Detection (IHB)/ Atrial fibrillation (AFIB)
- Movement Error
- Pressure Rating Indicator
- Battery Life Indicator
- Cuff Fit Error



Hypertension
CANADA

Voice Guided

One Button
Operation



Unlock the Potential
Download the FREE
A&D Heart Track
Companion App



Technical Specifications

| | |
|-------------------------------------|--|
| Measurement Method | Oscillometric |
| Memory | 90 measurements |
| Cuff Size | 22-42 cm (8.6-16.5 in) |
| Measurement Range | Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats/min |
| Measurement Accuracy | Pressure: ±3 mmHg, Pulse: ±5% |
| Power Supply | 4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TM-233) |
| Battery Life | Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries) |
| Monitor Dimensions | 14 x 11.8 x 6.4 cm (5.5 x 4.6 x 2.5 in) |
| Monitor Weight | 320 g (11.3 oz) without batteries |
| Display Size (WxH) | 6.2 x 7.2 cm (2.4 x 2.8 in) |
| Clinical Test | ISO81060-2: 2020 |
| Operating Conditions | 10°C to 40°C (+50°F to +104°F) / 15% to 85% RH / 800 to 1060 hPa |
| Storage/Transport Conditions | -20°C to 60°C (-4°F to +140°F) / 10% to 95% RH / 700 to 1060 hPa |

Packaging and Shipping Specifications

| | Dimensions (WxLxH) | Weight | Case Cube | UPC / GTIN |
|----------------------------|---|------------------|--|--------------------|
| Individual Unit Box | 13 x 10.5 x 19.5 cm (5.1 x 4.1 x 7.7 in) | 770 g (1.7 lb) | 0.003 m ³ (0.09 ft ³) | 0-93764-60457-7 |
| Case (10 units) | 28 x 55 x 22 cm (11 x 21.7 x 8.7 in) | 8.4 kg (18.5 lb) | 0.03 m ³ (1.2 ft ³) | 1-00-93764-60457-4 |

Optional Accessories

| Name | Order Code |
|------------------|------------|
| AC Power Adapter | TB-233 |

Compatible Replacement Cuffs- Packaged in individual box

| | |
|--|---------|
| Small Cuff 16-24 cm (6.3-9.4 in) | UA-289A |
| Medium Cuff 23-37 cm (9.0-14.6 in) | UA-290A |
| Large Cuff 31-45 cm (12.2-17.7 in) | UA-291A |
| Wide Range Cuff 22-42 cm (8.6-16.5 in) | UA-420A |
| Cuffs are not made with natural rubber latex | |

What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff 22-42 cm (8.6-16.5 in)
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 4 AA batteries