

A traditional blood pressure monitor that delivers quick and accurate readings at the touch of a button.

Product Benefits

- Easily track blood pressure over time with 60 measurement memory eliminating the need for written log.
- Displays the average of the last three measurements providing more consistent blood pressure data.

Comfort and Convenience

- One button operation for quick and easy functionality.
- Advanced inflation and deflation technology provides faster and more comfortable readings.
- Includes a carrying case and batteries for immediate use and portability.

The A&D Promise

- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting.
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value.



60 Measurement Memory

Pressure Rating Indicator

Battery Life Indicator



Hypertension CANADA

Irregular Heartbeat Detection (IHB)/ Atrial fibrillation (AFIB)

One Button Operation



Unlock the Potential
Download the FREE
A&D Heart Track
Companion App
for manual entry



Technical Specifications

Measurement Method	Oscillometric
Memory	60 measurements
Cuff Size	22-42 cm (8.6-16.5 in)
Measurement Range	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats/min
Measurement Accuracy	Pressure: ±3 mmHg, Pulse: ±5%
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	9.6 x 13.8 x 6.7 cm (3.8 x 5.4 x 2.6 in)
Monitor Weight	9.2 oz (260g) without batteries
Display Size (WxH)	3.5 x 5.3 cm (1.4 x 2.1 in)
Clinical Test	ISO81060-2: 2020
Operating Conditions	10°C to 40°C (+50°F to +104°F) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-20°C to 60°C (-4°F to +140°F) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Case Cube	UPC / GTIN
Individual Unit Box	12.5 x 9 x 18 cm (4.9 x 3.5 x 7.1 in)	620 g (1.4 lb)	0.002 m ³ (0.07 ft ³)	0-93764-60449-2
Case (10 units)	27 x 47.5 x 20 cm (10.6 x 18.7 x 7.9 in)	6.8 kg (14.9 lb)	0.03 m ³ (0.91 ft ³)	1-00-93764-60449-9

Optional Accessories

Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 16-24 cm (6.3-9.4")	UA-289A
Medium Cuff 23-37 cm (9.0-14.6")	UA-290A
Large Cuff 31-45 cm (12.2-17.7")	UA-291A
Wide Range Cuff 22-42 cm (8.6-16.5")	UA-420A
Cuffs are not made with natural rubber latex	

What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff 22-42 cm (8.6-16.5 in)
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 4 AA batteries