

Advanced, multi user Bluetooth® blood pressure monitor with enhanced features that provide trusted accuracy and reliable measurements every time.

Product Benefits

- Experience accurate readings with advanced cuff fit detection, avoiding errors from incorrect cuff placement.
- Ensure precise measurements with integrated movement detection, minimizing the risk of inaccurate results due to movement.
- Automatically records and tracks measurements with the highly rated Heart Track app to enable trending and sharing for better health management.
- Displays the average of the last three measurements providing more consistent blood pressure data.

Comfort and Convenience

- Up to 4 users provides convenience, cost savings, and enhanced data management capabilities for households with multiple individuals needing to monitor their blood pressure regularly.
- Advanced inflation and deflation technology provides faster and more comfortable readings.

The A&D Promise

- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting.
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value.



Advanced Averaging


60 Measurement Memory Per User

Movement Error

Pressure Rating Indicator

Battery Life Indicator

Multi User



Irregular Heartbeat Detection (IHB)/ Atrial fibrillation (AFIB)

Cuff Fit Error

One Button Operation



Hypertension CANADA



Unlock the Potential
Download the FREE
A&D Heart Track
Companion App



Technical Specifications

Measurement Method	Oscillometric
Memory	60 measurements per user
Cuff Size	22-42 cm (8.6-16.5 in)
Measurement Range	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats/min
Measurement Accuracy	Pressure: ±3 mmHg, Pulse: ±5%
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-233)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	14 x 11.8 x 6.4 cm (5.5 x 4.6 x 2.5 in)
Monitor Weight	10.4 oz (295g) without batteries
Display Size (WxH)	4.3 x 4.9 cm (1.7 x 1.9 in)
Clinical Test	ISO81060-2: 2020
Operating Conditions	10°C to 40°C (+50°F to +104°F) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-20°C to 60°C (-4°F to +140°F) / 10% to 95% RH / 700 to 1060 hPa -20°C to 60°Cv

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Case Cube	UPC / GTIN
Individual Unit Box	13 x 10.5 x 19.5 cm (5.1 x 4.1 x 7.7 in)	750 g (1.7 lb)	0.003 m ³ (0.09 ft ³)	0-93764-60450-8
Case (10 units)	28 x 55 x 22 cm (11 x 21.7 x 8.7 in)	8.2 kg (18.1 lb)	0.03 m ³ (1.2 ft ³)	1-00-93764-60450 5

Optional Accessories

Name	Order Code
AC Power Adapter	TB-233

Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 6-24 cm (6.3-9.4 in)	UA-289A
Medium Cuff 23-37 cm (9.0-14.6 in)	UA-290A
Large Cuff 31-45 cm (12.2-17.7 in)	UA-291A
Wide Range Cuff 22-42 cm (8.6-16.5 in)	UA-420A
Cuffs are not made with natural rubber latex	

What is In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff 22-42 cm (8.6-16.5 in)
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 4 AA batteries