

Advanced Multi User Bluetooth® blood pressure monitor with a suite of cutting-edge features that deliver Clinical Grade Accuracy and reliable measurements every time

Product Benefits

- Experience accurate readings with advanced cuff fit detection, avoiding errors from incorrect cuff placement
- Ensure precise measurements with integrated movement detection, minimizing the risk of inaccurate results due to movement
- TriCheck™ minimizes the impact of individual measurement fluctuations by automating 3 sequential measurements and averaging in accordance with clinical guidelines
- For additional peace of mind, device will check for cuff leaks to ensure accuracy
- Automatically records and tracks 90 measurements on device and unlimited on the highly rated Heart Track app to enable trending and sharing for better health management
- Displays the average of AM and PM measurements, in accordance with healthcare guidelines, to provide a more comprehensive assessment of fluctuations



Comfort and Convenience

- Up to 2 users provides convenience, cost savings, and enhanced data management capabilities for households with multiple individuals needing to monitor their blood pressure regularly
- Ability to adjust the resting time between measurements ensures more accurate and consistent readings
- Large, backlit display provides clear, easy-to-read measurements

The A&D Promise

- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value

Measurement Memory

Movement

Cuff Fit Error

TriCheck

Pressure Rating Indicator

Battery Life Indicator

User



Clock



Set Button

User/ Self Check Button

Start Button



Unlock the Potential
Download the FREE
A&D Heart Track
Companion App



Technical Specifications

Measurement Method	Oscillometric
Memory	90 measurements per user
Cuff Size	22-42 cm (8.6-16.5")
Measurement Range	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats / minute
Measurement Accuracy	Pressure: ± 3 mmHg, Pulse: $\pm 5\%$
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-304)
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 150 times R6P (manganese batteries)
Monitor Dimensions	Approx. 5.5 x 2.5 x 4.3" (13.4 x 6.3 x 10.9 cm)
Monitor Weight	Approx. 12.4 oz. (351g), excluding the batteries
Clinical Test	ISO81060-2: 2020
Operating Conditions	+50° F to +104° F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4 to +140° F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Individual Item Cube	UPC / GTIN
UA-810BLE Individual Unit (Box)	5.1 x 4.1 x 7.7 in (13 x 10.5 x 19.5 cm)	1.7 lb (750 g)	0.09 cu ft (0.003 cu meter)	0 93764 60452 2
UA-810BLE Case (10 units)	11 x 21.7 x 8.7 in (28 x 55 x 22 cm)	18.2 lb (8.3 kg)	1.2 cu ft (0.03 cu meter)	1 00 93764 60452 9

Optional Accessories

Name	Order Code
AC Power Adapter	TB:304
Cuff Holder	AX-CUF-HLD

Included In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff (8.6-16.5")
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 1 AC Adapter

Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 6.3-9.4" (16-24 cm)	UA-289A
Medium Cuff 9.0-14.6" (23-37 cm)	UA-290A
Large Cuff 12.2-17.7" (31-45 cm)	UA-291A
Wide Range Cuff 8.6-16.5" (22-42 cm)	UA-420A
Preformed Wide Range Cuff 8.6-16.5" (22-42 cm)	UA-372A
Cuffs are not made with natural rubber latex	



4622 Runway Boulevard
Ann Arbor, MI 48108
1-888-726-9965
andmedical.com