

Next Generation Multi User Bluetooth® blood pressure monitor with an enhanced suite of cutting-edge features, delivering unparalleled Clinical Grade Accuracy and consistently reliable measurements.

**Product Benefits**

- Experience accurate readings with advanced cuff fit detection, avoiding errors from incorrect cuff placement
- Ensure precise measurements with integrated movement detection, minimizing the risk of inaccurate results due to movement
- Device displays ambient room temperature to help end users understand the impact of temperature variations on their blood pressure
- TriCheck™ minimizes the impact of individual measurement fluctuations by automating 3 sequential measurements and averaging in accordance with clinical guidelines
- For additional peace of mind, device will check for cuff leaks and has dual sensors for calibration to ensure greater accuracy and reliability.
- Automatically records and tracks 120 measurements per user on device and unlimited on the highly rated Heart Track app to enable trending and sharing for better health management
- Displays the average of AM and PM measurements, in accordance with healthcare guidelines, to provide a more comprehensive assessment of fluctuations




**Comfort and Convenience**

- Up to 4 users provides convenience, cost savings, and enhanced data management capabilities for households with multiple individuals needing to monitor their blood pressure regularly
- Large, backlit display provides clear, easy-to-read measurements

**The A&D Promise**

- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value



Unlock the Potential  
Download the FREE  
**A&D Heart Track**  
Companion App



## Technical Specifications

<b>Measurement Method</b>	Oscillometric
<b>Memory</b>	120 measurements per user
<b>Cuff Size</b>	22-42 cm (8.6-16.5")
<b>Measurement Range</b>	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats / minute
<b>Measurement Accuracy</b>	Pressure: $\pm 3$ mmHg, Pulse: $\pm 5\%$
<b>Power Supply</b>	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-304)
<b>Battery Life</b>	Approx. 700 times LR6 (alkaline batteries) Approx. 150 R6P (manganese batteries)
<b>Monitor Dimensions</b>	Approx. 5.5 x 2.5 x 4.3" (13.9 x 6.4 x 11 cm)
<b>Monitor Weight</b>	Approx. 12.4 oz.(351 g), excluding the batteries
<b>Clinical Test</b>	ISO81060-2: 2020
<b>Operating Conditions</b>	+50° F to +104° F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
<b>Storage/Transport Conditions</b>	-4 to +140° F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa

## Packaging and Shipping Specifications

	<b>Dimensions (WxLxH)</b>	<b>Weight</b>	<b>Individual Item Cube</b>	<b>UPC / GTIN</b>
<b>UA-910BLE Unit</b>	6.5 x 3.5 x 8.3 in (16.6 x 9 x 21 cm)	1.8 lb (810 g)	0.11 cu ft (0.003 cu meter)	0 93764 60453 9
<b>UA-910BLE Case (10)</b>	13.9 x 19.9 x 8.7 in (35.4 x 50.6 x 22 cm)	19.7 lb (9 kg)	1.4 cu ft (0.04 cu meter)	1 00 93764 60453 6

## Optional Accessories

Name	Order Code
AC Power Adapter	TB-304
Cuff Holder	AX-CUF-HLD

## Included In The Box

- 1 Blood pressure monitor
- 1 Pre-formed Blood pressure cuff (8.6-16.5")
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 1 AC Adapter

### Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 6.3-9.4" (16-24 cm)	UA-289A
Medium Cuff 9.0-14.6" (23-37 cm)	UA-290A
Large Cuff 12.2-17.7" (31-45 cm)	UA-291A
Wide Range Cuff 8.6-16.5" (22-42 cm)	UA-420A
Preformed Wide Range Cuff 8.6-16.5" (22-42 cm)	UA-372A
Cuffs are not made with natural rubber latex	



4622 Runway Boulevard  
Ann Arbor, MI 48108  
1-888-726-9965  
andmedical.com