

**Advanced Multi User Bluetooth® blood pressure monitor with a suite of cutting-edge features that deliver Clinical Grade Accuracy and reliable measurements every time**

### Product Benefits

- Experience accurate readings with advanced cuff fit detection, avoiding errors from incorrect cuff placement
- Ensure precise measurements with integrated movement detection, minimizing the risk of inaccurate results due to movement
- TriCheck™ minimizes the impact of individual measurement fluctuations by automating 3 sequential measurements and averaging in accordance with clinical guidelines
- For additional peace of mind, device will check for cuff leaks to ensure accuracy
- Automatically records and tracks 90 measurements on device and unlimited on the highly rated Heart Track app to enable trending and sharing for better health management
- Displays the average of AM and PM measurements, in accordance with healthcare guidelines, to provide a more comprehensive assessment of fluctuations

### Comfort and Convenience

- Up to 2 users provides convenience, cost savings, and enhanced data management capabilities for households with multiple individuals needing to monitor their blood pressure regularly
- Ability to adjust the resting time between measurements ensures more accurate and consistent readings
- Large, backlit display provides clear, easy-to-read measurements



### The A&D Promise

- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value



Hypertension  
CANADA



**Unlock the Potential**  
Download the FREE  
**A&D Heart Track**  
Companion App  
for manual entry



## Technical Specifications

<b>Measurement Method</b>	Oscillometric
<b>Memory</b>	90 measurements per user
<b>Cuff Size</b>	22-42 cm (8.6-16.5")
<b>Measurement Range</b>	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats / minute
<b>Measurement Accuracy</b>	Pressure: $\pm 3$ mmHg, Pulse: $\pm 5\%$
<b>Power Supply</b>	4 x 1.5V batteries (R6P, LR6 or AA) or AC adapter (TB-304)
<b>Battery Life</b>	Approx. 700 times LR6 (alkaline batteries) Approx. 150 times R6P (manganese batteries)
<b>Monitor Dimensions</b>	Approx. 13.4 x 6.3 x 10.9 cm (5.5 x 2.5 x 4.3")
<b>Monitor Weight</b>	Approx. 12.4 oz. (351g), excluding the batteries
<b>Clinical Test</b>	ISO81060-2: 2020
<b>Operating Conditions</b>	10°C to 40°C ( +50° F to +104° F) / 15% to 85% RH / 800 to 1060 hPa
<b>Storage/Transport Conditions</b>	-20°C to 60°C ( -4° F to +140° F) / 10% to 95% RH / 700 to 1060 hPa

## Packaging and Shipping Specifications

	<b>Dimensions (WxLxH)</b>	<b>Weight</b>	<b>Individual Item Cube</b>	<b>UPC / GTIN</b>
<b>UA-810CNBLE Individual Unit (Box)</b>	13 x 10.5 x 19.5 cm (5.1 x 4.1 x 7.7 in)	750 g (1.7 lb)	0.003 m <sup>3</sup> (0.09 ft <sup>3</sup> )	0 93764 60458 4
<b>UA-810CNBLE Case (10 units)</b>	28 x 55 x 22 cm (11 x 21.7 x 8.7 in)	8.3 kg (18.2 lb)	0.03 m <sup>3</sup> (1.2 ft <sup>3</sup> )	1 00 93764 60458 1

## Optional Accessories

Name	Order Code
AC Power Adapter	TB:304
Cuff Holder	AX-CUF-HLD

### Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 16-24 cm (6.3-9.4" )	UA-289A
Medium Cuff 23-37 cm (9.0-14.6")	UA-290A
Large Cuff 31-45 cm (12.2-17.7")	UA-291A
Wide Range Cuff 22-42 cm (8.6-16.5")	UA-420A
Preformed Wide Range Cuff 22-42 cm (8.6-16.5")	UA-372A
Cuffs are not made with natural rubber latex	

## Included In The Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff (22-42 cm)
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 1 AC Adapter