

Advanced Multi User Bluetooth® blood pressure monitor with a suite of cutting-edge features that deliver Clinical Grade Accuracy and reliable measurements every time.



Hypertension
CANADA



Clinical Grade Accuracy

Features

- Up to 2 users - 90 memory
- Measurement Accuracy – Movement and Cuff Fit error detection
- Precision Check - Checks for cuff leaks to ensure accuracy
- Automatically records and tracks on the A&D Heart Track App
- Advanced Averaging – Last 3 measurements with rest periods
- Irregular Heartbeat / Atrial Fibrillation (AFib)
- Pressure Rating Indicator



Unmatched Reliability

With an industry-leading 5-Year warranty on monitor, we ensure consistent, trustworthy devices and measurements



Customer Satisfaction

Our dedicated service team offers seamless local assistance in English and French



Trusted Experts

Developed by A&D Medical, a pioneer in health monitoring technology

Technical Specifications

Measurement Method	Oscillometric
Memory	90 measurements per user
Cuff Size	22-42 cm (8.6-16.5")
Measurement Range	Pressure: 0-299 mmHg; Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg; Pulse: 40-180 beats/min
Measurement Accuracy	Pressure: ± 3 mmHg, Pulse: $\pm 5\%$
Power Supply	4 x 1.5V batteries (R6P, LR6 or AA), AC adapter (TB-304)
Battery Life	Approx. 700 times LR6 (alkaline batteries); Approx. 150 times R6P (manganese batteries)
Monitor Dimensions	(WxHxD) Approx. 13.9x6.4x11.2 cm (5.5x2.5x4.4 in)
Monitor Weight	Approx. 12.4 oz.(351 g), excluding the batteries
Wireless Communication	Bluetooth®: Ver.5.1LE BLP
Clinical Test	ISO81060-2: 2020
Operating Conditions Storage/Transport Conditions	10°C to 40°C (+50° F to +104° F) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-20°C to 60°C (-4 to +140° F) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Case Cube	UPC / GTIN
Individual Unit	13x10.5x19.5 cm (5.1x4.1x7.7 in)	750 g (1.7 lb)	0.003 m³ (0.09 ft³)	0 93764 60458 4
Case (x units)	28x55x22 cm (11x21.7x8.7 in)	8.3 kg (18.2 lb)	0.03 m³ (1.2 ft³)	1 00 93764 60458 1

Optional Accessories

Part Number	Description
UA-289A	Small Cuff 6.3-9.4 in (16-24 cm)
UA-290A	Medium Cuff 9.0-14.6 in (23-37 cm)
UA-291A	Large Cuff 12.2-17.7 in (31-45 cm)
UA-420A	Wide Range Cuff 8.6-16.5 in (22-42 cm)
UA-372A	Preformed Wide Range Cuff 8.6-16.5 in (22-42 cm)
TB-304	AC Adapter

What is in the Box

- 1 Blood pressure monitor
- 1 Blood pressure cuff (8.6-16.5 in)
- 1 Instruction manual
- 1 Quick start guide
- 1 Carrying case
- 1 AC Adapter



6695 Millcreek Drive, Unit 6
Mississauga, ON L5N 5R8, Canada
andmedical.com



Unlock the Potential
Download the **FREE**
A&D Heart Track
Companion App



All trademarks are property of their rightful owners. 09/2025.