



	<i>Easy-on PC</i> (2700-3)	<i>EasyOne Air</i> (2500-2A)
<b>Spirometer type</b>	PC spirometer	Portable & PC spirometer
<b>Test types</b>	FVC, FVL, SVC, MVV, Provocation	FVC, FVL, Tidal FVC, Tidal FVL, SVC & MVV
<b>EMR Connectivity</b>	Connectivity utilizing EasyOne Connect software for HL7 and XML interface	Connectivity utilizing EasyOne Connect software for HL7 and XML interface
<b>Spirometry incentives</b>	Real-time curves & graphs with EasyOne Connect software	Real-time curves on device & incentives available via Bluetooth to EasyOne Connect software
<b>Software</b>	EasyOne Connect software	EasyOne Connect software
<b>Cybersecurity</b>	EasyOne Connect – customer-managed encrypted databases	EasyOne Connect - customer-managed encrypted databases
<b>Technology</b>	Ultrasound TrueFlow Technology	Ultrasound TrueFlow Technology
<b>Calibration free</b>	Calibration not required, device maintains stability for a lifetime of use	Calibration not required, device maintains stability for a lifetime of use
<b>Power source</b>	USB connection	Rechargeable lithium-ion battery
<b>2019 ATS/ERS standards</b>	Meets all current guidelines	Meets all current guidelines
<b>Standalone testing</b>	Fully-integrated PC-connected spirometer (laptop/computer not included)	Portable or PC mode for multi-use application
<b>Mouthpiece</b>	nnd Spirette	nnd FlowTube
<b>Warranty</b>	2 or 3 years (depending on location)	2 or 3 years (depending on location)