

A simple and accurate blood pressure monitor providing measurements at the touch of a button and includes a large cuff for reliable measurements for people with larger arms.

Product Benefits

- Displays the average of all measurements providing more consistent blood pressure data
- Large cuff fits arms (12.2-17.7" / 31-45cm)
- Compatible with Small to Large cuffs (not included)
- Easily track blood pressure over time with 60 measurement memory eliminating the need for written log

Comfort and Convenience

- One button operation for quick and easy functionality
- Advanced inflation and deflation technology provides faster and more comfortable readings
- Adjustable pressure settings to ensure comfort
- Includes AC adapter for immediate use

The A&D Promise


- Clinical Grade Accuracy delivers measurements that are as precise as those taken in a medical setting
- Extended 5-year warranty on monitor provides long-term peace of mind and greater value



Measurement Memory

Pressure Rating Indicator

Battery Life Indicator



Unlock the Potential
Download the FREE
A&D Heart Track
Companion App
for manual entry



Irregular Heartbeat
Symbol (IHB)/ AFib

One Button
Operation

Technical Specifications

Measurement Method	Oscillometric
Memory	60 measurements
Cuff Size	31-45cm (12.2-17.7")
Measurement Range	Pressure: 0-299 mmHg, Systolic pressure: 60-279 mmHg Diastolic pressure: 40-200 mmHg, Pulse: 40-180 beats/min
Measurement Accuracy	Pressure: ±3 mmHg, Pulse: ±5%
Power Supply	AC Adapter (TB:233) included or 4 x 1.5V batteries (R6P, LR6 or AA) sold separately
Battery Life	Approx. 700 times LR6 (alkaline batteries) Approx. 200 times R6P (manganese batteries)
Monitor Dimensions	3.8 x 5.1 x 2.7" (9.6 x 13 x 6.8 cm)
Monitor Weight	8.5 oz (240 g) without batteries
Display Size	1.3 x 2.1" (3.4 x 5.3 cm)
Clinical Test	ISO81060-2:2013
Operating Conditions	+50°F to +104°F (10°C to 40°C) / 15% to 85% RH / 800 to 1060 hPa
Storage/Transport Conditions	-4°F to +140°F (-20°C to 60°C) / 10% to 95% RH / 700 to 1060 hPa

Packaging and Shipping Specifications

	Dimensions (WxLxH)	Weight	Case Cube	UPC / GTIN
Individual Unit Box	12.5 x 9 x 18 cm (4.9 x 3.5 x 7.1")	0.7 kg (1.6 lb)	0.002 m ³ (0.07 ft ³)	0-93764-60441-6
Case (10 units)	47.5 x 26.7 x 19.4 cm (18.7 x 10.5 x 7.6")	6.5 kg (14.3 lb)	0.03 m ³ (0.87ft ³)	1 00 93764 60441 3

Optional Accessories

Name	Order Code
AC Power Adapter	TB-233

Compatible Replacement Cuffs- Packaged in individual box

Small Cuff 16-24 cm (6.3-9.4")	UA-289A
Medium Cuff 23-37 cm (9.0-14.6")	UA-290A
Large Cuff 31-45 cm (12.2-17.7")	UA-291A
Wide Range Cuff 22-42 cm (8.6-16.5")	UA-420A
Cuffs are not made with natural rubber latex	

What is In The Box

- 1 Monitor
- 1 Blood pressure cuff Large cuff fits arms 31-45cm (12.2-17.7")
- 1 Instruction Manual
- 1 Quick Start Guide
- 1 AC Adapter

Health Canada License - 92815



6695 Millcreek Drive, Unit 6
Mississauga, ON L5N 5R8, Canada
If you want to place your order, please email us to orders@autocontrol.com or call our customer service team
1-800-461-0991
andmedical.com